International assembly and repair
instructions

Xlevator

Languages:

EN, DE, SL

fEX Engineering GmbH

Höltenweg 49
48155 Münster
Germany

Third edition

Publication date 22/07/2024

Foreword

Thank you for purchasing an **Xlevator** chair bracket from fEX Engineering GmbH.

The Xlevator guarantees that the powered **Xpedition** stair chair made by Stryker is properly secured. Power-assisted linear movement ensures ergonomic handling when removing or stowing it in a rescue vehicle.

WARNING

Improper handling or any unintended use may cause incur risk or cause damage. Please therefore read these instructions for use carefully and follow them closely. Retain them for reference at all times. To prevent personal injury or damage to property, please observe all safety precautions.

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster
Germany
engineering@fex.group
+49 (0)251 508 537 -13

Contents

1	How to use these Assembly and repair instructions	1
1.1	Feedback on these Instructions for use	1
1.2	Explanation of the warnings used in these Instructions for use	1
2	Product description	1
2.1	Intended use.....	1
2.2	Product overview	2
2.3	How it works.....	3
2.4	Symbols used.....	3
2.5	Use in combination with other products	4
2.6	Delivery.....	4
2.7	Modifications.....	4
2.8	Requirements for storage, handling and environmental conditions	4
2.9	Product service life	4
2.10	Warranty.....	4
3	Installation on a rear door or wall	5
3.1	Requirements for personnel involved	5
3.2	General	5
3.3	Screw connections.....	6
3.4	Gas spring and oil brake	7
3.5	Wheel mount.....	8
3.6	To check for correct assembly.....	8
4	Maintenance.....	9
4.1	Maintenance.....	10
4.2	Maintenance by technicians.....	10
4.3	Gas spring, oil brake and cylinder connector	10
4.4	Replacing the gas spring.....	11
4.5	Replacing the oil brake	12
4.6	Replacing the locking unit	13
4.7	Lubricating the guide rails	14
5	Dismantling, recycling, disposal	15
5.1	Dismantling.....	15
5.2	Recycling, disposal.....	15
5.3	Disposal through fEX Engineering GmbH	15
6	List of individual parts	16

1 How to use these Assembly and repair instructions

If you have any questions about the contents of these instructions for use or how to use the product, please contact us. You can request a paper copy of these instructions for use or download the current, digital version from our website:

<https://www.fex.group/downloads>

Should you need further assistance after carefully reading these Instructions for use and additional information, please contact your authorised distributor or contact us directly.

1.1 Feedback on these Instructions for use

Your opinion is important to us. Please let us know if there is anything you would change about these Instructions for use (engineering@fex.group). We will take your feedback into account in the next version of these Instructions for use.

1.2 Explanation of the warnings used in these Instructions for use

Signal word	Description
 WARNING	Potentially hazardous situation resulting in severe injury or death.
 CAUTION	Potential hazard resulting in minor or moderate injury.
NOTE	User error can damage the product.

Table 1: Explanation of the warning notices used in these Instructions for use

2 Product description

2.1 Intended use

The Xlevator is designed as a professional securing device for the powered **Xpedition** stair chair made by *Stryker*. The system is intended to raise and lower the *Xpedition Chair* between its stowed position and ground level. Power-assisted linear movement ensures ergonomic handling when removing or stowing it in a rescue vehicle.

2.2 Product overview


Table 2: Product overview

2.3 How it works

The Xlevator is an accessory for the powered *Xpedition* stair chair made by *Stryker*. It consists of a frame, a wall plate and movable carriage and a chair plate. These two elements are connected to one another with a linear movement, and can therefore be moved vertically. A gas spring and an oil brake are installed between the chair plate and the wall plate. The gas spring supports upward movement so that the user need not exert a lot of strength to stow the *Xpedition* in the emergency services vehicle. In contrast, the oil brake ensures that the chair plate rises in a controlled manner on empty trips. The built-in locking mechanism enables the *Xpedition* to be secured in incremental vertical positions. The handle bar of the *Xpedition* locks behind the bar clasp at the top of the chair plate. The wheels fit into two wheel mounts at the bottom of the chair plate. The *Xpedition* is then secured with two straps.

2.4 Symbols used

Figure 1: Labels

Figure 2: Caution: Crush hazard

Symbol	Meaning
	Manufacturer's logo
	Manufacturer, address and contact details
	Conformity marking
	Read Instructions for use
	Caution: Crush hazard
S/N	2024 03 06 001 Year of construction/Month/Date/Batch number

Table 3: Symbols used and their meanings

2.5 Use in combination with other products

The Xlevator is designed exclusively for carrying and securing the powered Xpedition stair chair made by *Stryker*. Only this product can be stowed and transported with adequate safety. The screw connections in the vehicle must withstand a force of 5000 N.

2.6 Delivery

The product is delivered in non-sterile, non-disinfected condition.

An electronic copy of the Instructions for use can be found via the following link:

<https://www.fex.group/downloads/>

A sticker on the system shows a QR code and the above link, both of which will take you to the electronic copy of the Instructions for use. This sticker must be removed before use.

2.7 Modifications

Unauthorised modifications to this product are prohibited.

2.8 Requirements for storage, handling and environmental conditions

Category	Information
Permitted humidity limits	5% - 98%
Permitted temperature limits (vehicle & storage)	10°C to +50°C
Maximum storage time before use. Cylinders can remain in the same position for a maximum of 6 months.	6 months

Table 4: Storage during transport, before and after use

2.9 Product service life

Under normal usage conditions and with appropriate maintenance, the service life is **seven years**.

The date for the product's service life is based on the delivery date from fEX Engineering.

2.10 Warranty

The Xlevator is subject to a statutory warranty of **two years** under EU and German law.

The oil brake, gas spring and straps are excluded from the above warranty. A **one-year** warranty applies to these components.

3 Installation on a rear door or wall

3.1 Requirements for personnel involved

Installation, maintenance and repair can be carried out by the operator and his personnel and by vehicle manufacturers if the instructions below are followed. It is the responsibility of the technical manager to appoint a qualified person to do the work.

3.2 General

The Xlevator is attached to the vehicle using screw connections. For this purpose, the screw connections shown in Figure 3 must be present. These screw connections must be able to withstand a total tensile force of 5000 N.


Figure 3: Screw connections

⚠ CAUTION

A torque wrench must be used for assembly and the specified torques must be observed.

⚠ WARNING

Medium-strength screw locking is strongly recommended. The required curing time depends on the information given by the screw locking manufacturer.

The Xlevator is delivered in non-disinfected condition. You must disinfect the system before first use. Instructions for disinfection can be found in the separate document "*Instructions for use - Xlevator*".

3.3 Screw connections

Figure 5: Positioning the wall plate


Figure 4: Position M8 screws


Figure 6: Mounting the brackets

Follow the steps below:

1. Press the locking button (1) and pull the wall plate (2) upwards so that the through holes (3/4) line up with the riveting nuts (Fig. 4) (the system may also need to be lined).
2. Apply a drop of medium-strength screw locking to the first (lowest) threads on **all** screws.
3. Fasten the wall plate (2) to the vehicle using the four M8x16 (round head) screws (4) and a tightening torque of 22 Nm (Fig. 5).
4. Place the bracket (5) on the wall plate (2) and fasten it with two M6x16 (round head) screws (6) each using a tightening torque of 13 Nm (Fig. 6).
5. Repeat this for the second bracket.

3.4 Gas spring and oil brake

⚠ CAUTION

The gas spring is pre-tensioned when it is delivered. Follow the steps in the instructions to avoid injury or damage.

Fig. 7: Position of cylinder


Fig. 8: Cable retainer

Fig. 9: Cylinder connector

Fig. 10: Assembly aids

Follow the steps below:

1. Press the locking button (1) and pull the chair plate upwards so that the through holes for the cylinders line up with the bracket through holes (2) (Fig. 7).
2. Cut the cable retainer (3) halfway up the system (Fig. 8).
3. Connect the cylinders and the brackets using the bolts (4) and the self-locking M6 nut (5) with a tightening torque of 11 Nm (Fig. 9).
4. Press the locking button (1) again, pull the chair plate down a few notches and remove the assembly aids (6) (Fig. 10).

3.5 Wheel mount

The Xlevator has one left and one right wheel mount (Fig. 10) (1). The wheel mounts are symmetrical. On the underside of the wheel mounts, one edge is right-angled and the other is rounded.


Figure 11: Left wheel mount


Figure 12: Screw connections for wheel mount

Follow the steps below:

1. Place the left (Fig. 11, 1) and the right wheel mount in the chair plate (Fig. 12).
2. The right-angled inner edge marked in red (Fig. 11) is fitted flush. The opposite outer edge is rounded to fit over the weld seam. The right wheel mount is a mirror image of the left.
3. Apply a drop of medium-strength screw locking to the first (lowest) threads on **all** screws.
4. Attach the new wheel mount (1) using the three new M8x12 countersunk head screws (2) with a tightening torque of 8 Nm.

3.6 To check for correct assembly

WARNING

Use the check list below to make sure the system is assembled or installed correctly. If any of these points do not apply, the product must be reassembled or reinstalled.

- All necessary screws must be used and no through holes should be visible: four M8x16 (flat head), four M6x16 (raised head), six M8x10 (round head)
- All screws should be flush and tightly screwed (use a torque wrench for assembly).
- Unhindered up and down movement
- The gas spring supports lifting
- The oil brake controls the speed when moving upwards (empty trip)
- The system can be locked into different positions
- The straps are undamaged and working

Following the proper assembly, the Xlevator is ready for immediate use.

4 Maintenance

The operator and his staff can maintain the product by following the instructions. Unauthorised repairs to this product are prohibited. The operator is permitted to replace certain components; these are divided into two groups:

1. Components that can be replaced by **users (laymen)** (instructions in the separate document ***Instructions for use – Xlevator***), e.g. straps
2. Components that must be replaced by **technicians** (Section 4.2 *Maintenance by technicians*), e.g. the gas spring

It is the operator's responsibility to employ suitable, qualified and trained personnel carry this out. A video on how to maintain the Xlevator is available on the manufacturer's homepage.

 WARNING

Medium-strength screw locking is strongly recommended. The required curing time depends on the information given by the screw locking manufacturer.

 WARNING

Only replacement parts made by fEX Engineering GmbH may be used. Using other parts may be dangerous. Repairs are prohibited.

The Instructions for use are decisive for all work on the product. Once damaged components have been replaced, the system can be put back into use.

fEX Engineering GmbH refuses to supply spare parts for, or carry out repairs on, any products that have reached the end of their service life.

4.1 Maintenance

The guide rails must be lubricated as required. This depends on various factors, including the frequency/intensity of their cleaning, but must take place **at least once a year** (Section 8.4.5 *Lubricating the guide rails*).

A visual inspection must be made every week (instructions in the document *Instructions for use – Xlevator*). Faulty components should be replaced as necessary.

4.2 Maintenance by technicians

The parts listed below can be replaced by trained technicians (training video) with original spare parts. These can be ordered upon request.

- Cylinder connector (Item no. fEX0140ZV)
- Gas spring (290 Nm) (Item no. fEX0140ZF)
- Oil brake (2,7 mm) (Item no. fEX0140OB)
- Locking unit (Item no. fEX0140VE)
- Locking button (yellow) (Item no. fEX0140KG)
- Locking button (red) (Item no. fEX0140KR)

WARNING

Only replacement parts made by fEX Engineering GmbH may be used. Using other parts may be dangerous.

4.3 Gas spring, oil brake and cylinder connector

Only a trained technician (training video) is permitted to replace these parts. Both the gas spring and the oil brake must be installed with the piston rod pointing downwards. The gas spring and the oil brake can be differentiated by the length of their cylinders. **The cylinder of the gas spring is approx. 6 cm longer than the cylinder of the oil brake (Fig. 13).**

The cylinder connector is clamped around both cylinders with a screw and must be positioned almost at the bottom.


Figure 13: Identifying the cylinder

4.4 Replacing the gas spring

Only a trained technician (training video) is permitted to replace the gas spring.

⚠ CAUTION

The gas spring is pre-tensioned when it is delivered. Follow the steps in the instructions to avoid injury or damage.

Follow the steps below:

1. Move the system to its lowest position and lock into place.
2. Loosen the screw of the cylinder connector (1) (Fig. 14) and pull it down so that it hangs freely over the piston rods.
3. Identify the gas spring (see Section 3.3 *Gas spring, oil brake and cylinder connector*)
4. Place **both** assembly aids (2) (Fig. 15) over the piston rod on the gas spring.
5. Let the system move up as far as possible (assembly aid blocks the last part of the travel path).
6. Loosen the nut (3) on the upper bracket (Fig. 16) and remove the bolt (4).
7. Move the system down one notch and let it snap into place.
8. Repeat step 6 on the lower bracket.
9. Remove the gas spring from below.
10. Insert the new gas spring with the piston rod facing downwards. Guide the cylinder connector over the piston rod.
11. Connect the cylinder and the **lower** bracket using the bolt and the self-locking M6 nut with a tightening torque of 11 Nm.
12. Press the locking button and pull the chair plate upwards so that the through holes of the cylinders line up with the bracket through holes (Fig. 16).
13. Connect the gas spring and the **upper** bracket (Fig. 16) with the bolt (4) and the self-locking M6 nut (3) using a tightening torque of 11 Nm.
14. Pull the cylinder connector over both cylinders (Fig. 14) and fasten it to the M5x10 screw (round head) (1) using a tightening torque of 5 Nm.
15. Use the product to check it is working properly.

Fig. 14: Cylinder connector


Fig. 15: Assembly aids

Fig. 16: Upper bracket

4.5 Replacing the oil brake

Only a trained technician (training video) is permitted to replace the oil brake.

Follow the steps below:

1. Move the system to the middle position and lock into place.
2. Loosen the screw of the cylinder connector (1) (Fig. 17) and pull it down so that the cylinder connector hangs freely over the piston rods.
3. Identify the gas spring (see Section 8.5.1 *Gas spring, oil brake and cylinder connector*).
4. Loosen the nut (2) on the upper bracket (Fig. 18) and remove the bolt (3).
5. Loosen the nut on the lower bracket and remove the bolt.
6. Remove the oil brake from below.
7. Insert the new oil brake with the piston rod facing downwards. Guide the cylinder connector over the piston rod.
8. Connect the cylinder and the **lower** bracket using the bolt and the self-locking M6 nut with a tightening torque of 11 Nm.
9. Pull the oil brake upwards by hand so that the through holes of the cylinder line up with the bracket through holes.
10. Connect the cylinder and the **upper** bracket (Fig. 18) with the bolt (3) and the self-locking nut (2) using a tightening torque of 11 Nm.
11. Pull the cylinder connector (Fig. 17) over both cylinders and fasten it with the M5x10 screw (round head) (1) using a tightening torque of 5 Nm.
16. Use the product to check it is working properly.

Fig. 17: Cylinder connector

Fig. 18: Upper bracket

4.6 Replacing the locking unit

Only a trained technician (training video) is permitted to replace the locking unit.

CAUTION

The chair plate must be in the highest position because the system moves upwards automatically. Follow the steps in the instructions to avoid injury or damage.

Follow the steps below:

1. Move the system to the middle position and let it lock into place.
2. Place **both** assembly aids (2) (Fig. 15) over the piston rod on the gas spring.
3. Press the locking button and let the system move up until the assembly aids fit.
4. Remove the handle by loosening the four M5x12 screws on the back of the chair plate.
5. Put your hand behind the locking unit and hold it firmly.
6. Remove the four countersunk head screws.
7. Turn the entire locking unit by approx. 45° and remove the unit.
8. Apply a drop of medium-strength screw locking to the new countersunk head screws.
9. Take the replacement unit and rotate it into position under the button. The button may need to be turned slightly so that the unit is flush with the plate.
10. Press the unit onto the plate and turn the unit back (approx. 45°).
11. Screw in the accompanying M6x16 (countersunk head screw) screws using a tightening torque of 10 Nm.
12. Move the system back to the middle position and let it lock into place.
13. Remove the two assembly aids.
14. Then check that the lock works in different positions.

4.7 Lubricating the guide rails

The rails should be lubricated as necessary, but **at least once a year**. As foreign bodies can get into the rails, these must be cleaned first. An NLGI Class 2 or similar lubricant must be used.

Follow the steps below:

1. Move the system to the lowest position.
2. Clean the rails with a dry cloth.
3. Lubricate both track surfaces on the rails (Fig. 19 and Fig. 20).
4. Repeat steps 2 and 3 for all four rails.
5. Move the system up and down to spread the lubricant.
6. Do **not** remove excess lubricant from inside the rails.


Figure 19: Wall plate rails


Figure 20: Chair plate rails

5 Dismantling, recycling, disposal

5.1 Dismantling

The product must be cleaned and disinfected before being dismantled. Components that have been dismantled or damaged may have sharp edges.

The operator and his personnel may dismantle the product. It is the operator's responsibility to have qualified personnel carry this out. The product is dismantled in the reverse order of its assembly.

5.2 Recycling, disposal

When disposing of the product, always follow **current local** recommendations, and/or statutory provisions on environmental protection and risk relating to recycling or the disposal of devices at the end of their useful life.

5.3 Disposal through fEX Engineering GmbH

If the product has to be taken out of service due to damage or other reasons, the Xlevator (and all other fEX Engineering GmbH products) can be sent to the company's address for free disposal.

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5.3 Odstranjevanje prek podjetja fEX Engineering GmbH