International Instructions for Use **Handle Extension**

Language:

EN, DA, DE,
ES, FI, FR, HU,
IT, NL, NO, SL,
SV

fEX Engineering GmbH

Höltenweg 49
48155 Münster
Germany

Fourth edition

Publication date 22/07/2024

Foreword

Thank you for purchasing a **handle extension** from fEX Engineering GmbH.

The handle extension is designed as an accessory for the *Power-PRO XT* wheeled stretcher made by *Stryker*. The handle bar is designed for to be in a raised position at the foot-end of the stretcher. It is specifically intended to help tall users by making it easier to manoeuvre when using the stretcher. The handle bar can also be removed to allow patients to be moved unimpeded, for example when treating the lower limbs or to use the maximum lying surface.

WARNING

Improper handling or any unintended use may cause incur risk or cause damage. Please therefore read and carefully follow these Instructions for use. Retain them for reference at all times. To prevent personal injury or damage to property, please observe all safety precautions.

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster
Germany
engineering@fex.group
+49 (0)251 508 537 -13

Contents

1	How to use these Instructions for use.....	1
1.1	Feedback on these Instructions for use	1
1.2	Explanation of the warning notices used in these Instructions for use	1
1.3	Abbreviations	2
2	Intended purpose	3
2.1	Intended use.....	3
2.2	Intended medical purpose	3
2.3	Clinical benefits	3
2.4	Indications	3
2.5	Contraindications	3
2.6	Intended target patient group	3
2.7	Intended use with other products.....	3
2.8	Intended users.....	4
2.9	Intended environment for use	4
3	Important safety notice.....	5
3.1	What to do in the event of damage	5
3.2	Reporting obligation.....	5
4	Product description	6
4.1	Product overview	6
4.2	Symbols used.....	7
4.3	How it works.....	7
4.4	Use in combination with other products	8
4.5	Delivery.....	8
4.6	Modifications.....	8
4.7	Requirements for storage, handling and environmental conditions	8
4.8	Product service life	8
5	Assembly and putting into service	9
5.1	Requirements for personnel involved.....	9
5.2	General	9
5.3	Assembly and installation.....	9
5.4	To check for correct assembly.....	10

6	Using the product.....	11
6.1	Operator's duty of care	11
6.2	User's duty of care.....	11
6.3	Instructions for safe operation	11
6.4	Quick-release fasteners	12
6.5	How to remove and insert the bar	13
7	Reprocessing: cleaning, disinfection	14
7.1	Cleaning the product	14
7.2	Disinfectants	14
7.3	Disinfecting the product.....	15
7.4	Final treatment.....	15
8	Maintenance.....	16
8.1	Maintenance.....	16
8.2	Inspection	16
8.3	Damage.....	16
8.4	Troubleshooting	17
8.5	Maintenance and repairs	18
9	Dismantling, recycling, disposal	19
9.1	Dismantling.....	19
9.2	Recycling, disposal.....	19
9.3	Disposal through fEX Engineering GmbH	19
10	Risk management	20
10.1	Foreseeable misuse	20
10.2	Implemented risk control measures	20
10.3	Residual risk.....	20
11	Performance data.....	21
11.1	Product conformity	21
12	Technical data.....	22
12.1	List of individual parts	22

1 How to use these Instructions for use

If you have any questions about the contents of these Instructions for use or how to use the product, please contact us. You can request a paper copy of these instructions for use or download the current, electronic version from our website:

<https://www.fex.group/downloads>

Should you need further assistance after carefully reading these Instructions for use and additional information, please contact your authorised distributor or contact us directly.

1.1 Feedback on these Instructions for use

Your opinion is important to us. Please let us know if there is anything you would change about these instructions for use (engineering@fex.group). We will analyse your feedback and take it into account in the next version of these Instructions for use.

1.2 Explanation of the warning notices used in these Instructions for use

Signal word	Description
 WARNING	Potentially hazardous situation resulting in severe injury or death.
 CAUTION	Potential hazard resulting in minor or moderate injury.
NOTE	User error can damage the product.

Table 1: Explanation of the warning notices used in these Instructions for use

1.3 Abbreviations

Abbreviation	Definition
BfArM	Bundesinstitut für Arzneimittel und Medizinprodukte (Federal Institute for Drugs and Medical Devices)
g	Gravitational constant
MDR	Medical Device Regulation
MPSV	Medizinprodukte-Sicherheitsplanverordnung (Medical Devices Safety Plan Ordinance)
Nm	Newton meter
PA	Polyamide
PL	Patient left
PR	Patient right
UDI	Unique Device Identification
UDI-DI	Unique Device Identification - Device Identifier
UNC	Unified coarse thread

Table 2: Abbreviations

2 Intended purpose

2.1 Intended use

The handle extension is only designed for pulling, pushing and steering, i.e. for manoeuvring the *Power-PRO XT* stretcher. Lifting the stretcher by the handle extension is prohibited.

2.2 Intended medical purpose

The handle extension is a Class I, mobile, non-invasive, non-active medical device.

2.3 Clinical benefits

The handle extension is classed as an accessory and has no inherent clinical benefit.

It is designed to help the *Power-PRO XT* wheeled stretcher made by *Stryker* to fulfil its intended purpose.

2.4 Indications

No indications have been specified for the handle extension.

2.5 Contraindications

There are no known contraindications for the handle extension.

2.6 Intended target patient group

The patient population is actually using the *Power-PRO XT* wheeled stretcher made by *Stryker*.

No direct contact with the patient is anticipated under normal circumstances.

2.7 Intended use with other products

The handle extension is intended exclusively as an accessory for the *Power-PRO XT* medical device made by *Stryker*.

2.8 Intended users

The handle extension is intended for use in a professional healthcare environment. The users, and the qualifications they require, are specified accordingly.

User groups	Nursing staff, emergency services personnel, doctors (primary users)	Medical technicians, technical managers, cleaning staff (secondary users)
Training, certification	Medically trained staff, e.g. emergency paramedics, emergency medical technicians, rescue workers, emergency services support staff, doctors, or qualified nurses	Medical technicians (training/academic qualifications) or equivalent, or accreditation as a cleaner in a professional healthcare environment

Table 3: Intended users

2.9 Intended environment for use

The handle extension is intended for use by primary users in the following environments:

- Use in proximity to patients and therefore in a professional healthcare environment
- Use in a mobile environment, e.g. in emergency services vehicles
- Use outdoors for rescue operations

The handle extension does not require any specific physical settings.

3 Important safety notice

NOTE

Please read and observe the instructions in these Instructions for use to use this system safely.

CAUTION

If a wheeled stretcher can no longer be used because of an accident or similar incident, the same applies to the handle extension fitted to it.

WARNING

The handle extension must be fully locked. Otherwise, there is a risk of the handle bar detaching from the stretcher, and the user losing control while manoeuvring the stretcher, or the handle bar becoming a projectile in the event of a traffic accident.

3.1 What to do in the event of damage

The handle extension must be withdrawn from use if it is damaged (Section 8.3 *Damage*). Repairing damaged components is prohibited – these must be replaced (Section 8.5 *Maintenance and repairs*).

Visible damage must be immediately reported to the technical manager.

3.2 Reporting obligation

All serious incidents which occur in connection with the product must be reported to the manufacturer and the competent authority in the EU Member State where you are established.

In Germany, the competent authority is BfArM. Up-to-date contact details can be found on the BfArM website: <https://www.bfarm.de>

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster
Germany
+49 (0)251 508 537 -13
engineering@fex.group

4 Product description

4.1 Product overview


Figure 1: Handle extension

- 1: Handle bar
- 2: Clamp assembly (PL)
- 3: Clamp assembly (PR)
- 4: Quick-release fastener

4.2 Symbols used

Figure 2: Notices on the product

Symbol	Meaning
	Manufacturer's logo
	Manufacturer, address and contact details
	Warning
	Medical device
	Conformity marking
	Read Instructions for use
	UDI - Unique Device Identification: (00) UDI; (10) Batch no.; (11) Date of manufacture

Table 4: Symbols on the product

4.3 How it works

The handle extension is an accessory for the *Power-PRO XT* wheeled stretcher made by *Stryker* and it helps this stretcher to fulfil its intended purpose.

This handle extension raises the handle bar which makes it easier to manoeuvre the stretcher, especially for tall users. It is also possible to remove the handle bar to reposition patients unimpeded. The handle bar can also be removed to maximise the lying surface or when treating the lower limbs. The handle bar is removed and locked via two lockable quick-release fasteners.

4.4 Use in combination with other products

The handle extension is designed exclusively for use with the *Power-PRO XT* wheeled stretcher made by *Stryker*. The product can be installed on this wheeled stretcher with positive locking.

The conditions for adequately safe use are only guaranteed when used in this way.

4.5 Delivery

The handle extension is delivered as a complete assembly, in non-sterile, non-disinfected condition.

An electronic copy of the Instructions for use can be found via the following link:

<https://www.fex.group/downloads/>

A sticker on the system shows a QR code and the above link, both of which will take you to the electronic version of the Instructions for use. This sticker must be removed before use.

4.6 Modifications

Unauthorised modifications to this product are prohibited.

4.7 Requirements for storage, handling and environmental conditions

The product requires no specific storage or handling conditions.

4.8 Product service life

Under normal usage conditions and with proper maintenance, the service life is **7 years**.

The date for the product's service life is based on the delivery date from fEX Engineering.

5 Assembly and putting into service

5.1 Requirements for personnel involved

The operator and his personnel can assemble the product by following the instructions below. It is the operator's responsibility to have qualified personnel carry this out.

5.2 General

The handle extension is supplied as a complete assembly, i.e. the inner clamp (PL and PR), outer clamp (PL and PR), handle bar and quick-release fasteners are packaged together. To ensure traceability, in accordance with the MDR, this assembly must be installed as it is delivered.

 CAUTION

A torque wrench must be used for assembly and the specified torques must be observed.

 WARNING

Medium-strength screw locking is strongly recommended. The required curing time depends on the information given by the screw locking manufacturer.

The handle extension is delivered in non-disinfected condition. You must disinfect the handle extension before first use. For disinfection instructions, see Section 7 *Reprocessing: cleaning and disinfection*.

5.3 Assembly and installation

The handle extension is designed exclusively for use with the *Power-PRO XT* wheeled stretcher made by *Stryker*. Positive locking and safe use is only guaranteed with this wheeled stretcher.

The system consists of a handle bar and four clamps. The clamps are screwed together in pairs on the patient left foot-end and the patient right foot-end (Figures 3 and 4). The two halves of each pair form a positive lock around the handle bars of the *Power-PRO XT*.


Figure 3: Clamp PR


Figure 4: Stryker assembly at foot-end

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Fit the *outside clamp PR* (1) from the outside onto the handle bar (6) (red area) of the assembly at the foot-end of the wheeled stretcher.
3. Fit the counterpart *inner clamp PR* (2) from the inside onto the handle bar (6) (red area).
4. Insert the four screws (3) UNC #10-24 x 1" and tighten them by hand.
5. Then tighten all four screws with a tightening torque of 7 Nm.
6. Repeat steps 2 and 5 as required for the other side.
7. Open the quick-release fasteners (4) and insert the handle bar (5) steadily.
8. Lock the quick-release fasteners (4).

5.4 To check for correct assembly

WARNING

Use the check list below to make sure the handle bar system is assembled or installed correctly. If any of these points do not apply, the product must be reassembled or reinstalled.

- There are four screws to attach on each side.
- All screws should be flush and tightly screwed (use a torque wrench for assembly).
- The gap between the inner clamp and the outer clamp must be completely closed by the sealing material.
- The handle bar cannot be removed without unlocking it.

Once assembled and installed properly, the handle extension is ready for immediate use

6 Using the product

6.1 Operator's duty of care

To ensure safe operation, the operator must meet the following obligations:

- Personnel must be trained to use the handle extension installed
- The Instructions for use must be read and followed
- Products must be assembled or installed properly
- Only suitable, qualified, trained specialist personnel may use the product
- All serious incidents (damage, injuries, etc.) must be reported to the manufacturer and the competent authorities

6.2 User's duty of care

To ensure safe operation, the user must meet the following obligations:

- The Instructions for use must be read and followed
- The product must be used as intended
- Damage must be reported to the technical manager

6.3 Instructions for safe operation

- Check the quick-release fastener is working
- Only use the product for its intended purpose
- Do not use to lift the stretcher
- Do not hang objects on the handle bar
- Do not use as a seat

6.4 Quick-release fasteners

The handle bar is held in place by fixed quick-release fasteners and can be unlocked manually by pulling up the locking pin (circled in red). To do this, pull the red button and turn it 30° to hold it in the locking position (6 locking positions). Pull and turn it again by 30° to drop the locking pin and lock the handle bar.

Figure 5: Quick-release fastener


Figure 6: In locked position

6 locking positions

Figure 7: Locking mechanism


Figure 8: Open position

6.5 How to remove and insert the bar

NOTE

The angle of the handle bar may cause it to get stuck when removing or reinserting it.

To remove the handle bar:

1. Open both quick-release fasteners. When the quick-release fasteners are open, the red button is clearly “open” (Section 6.4 *Quick-release fasteners*).
2. Hold the handle bar at both ends above the clamp (Figure 9, circled in red). Pull the bar out evenly on both sides.

To fit the handle bar:

1. Hold the bar at both ends above the clamp (Figure 9, circled in red).
2. Insert the bar evenly into the recesses with small movements.
3. Lock both quick-release fasteners. The quick-release fasteners are locked when the red buttons are clearly “down” (Section 6.4 *Quick-release fasteners*).

In both cases, be sure to move the handle bar evenly, or the handle bar may get stuck. If this happens, remove the bar again and reinsert it.


Figure 9: Clamp recess

⚠ WARNING

The handle extension must be fully locked. Otherwise, there is a risk of the handle bar detaching from the stretcher, and the user losing control while manoeuvring the stretcher, or the handle bar becoming a projectile in the event of a traffic accident.

7 Reprocessing: cleaning, disinfection

The handle extension is a reusable medical device. Cleaning and disinfection must follow the process in the Instructions for use, and be carried out after any possible soiling, or at least once a month.

We recommend cleaning and inspecting the handle extension at the same time. If possible, combine with existing weekly routines (e.g. when cleaning the vehicle).

7.1 Cleaning the product

The product can be cleaned using high pressure.

1. Remove the handle bar
2. Wipe the product's surfaces with a solution of mild soap and water.
3. Dry all product surfaces thoroughly.

7.2 Disinfectants

Use the following agents for disinfection:

- Quaternary cleaners (active ingredient: ammonium chloride) containing less than 3% glycol ether
- Chlorinated bleaching agents (5.25 % – less than 1 part bleaching agent to 100 parts of water)
- Phenolic cleaners (active ingredient: o-phenylphenol)
- ≤ 21% isopropanol

Disinfectants made from quaternary ammonium compounds must be used in accordance with the recommended instructions for use. Quaternary cleaners and disinfectants must be used with extreme care on plastic surfaces.

WARNING

Disinfectants containing any of the following basic ingredients are not approved as they can damage the system:

- Hydrogen peroxide
- (Per-/Chlor-) acetic acid

7.3 Disinfecting the product

1. Remove any visible contaminants from the product (Section 7.1 *Cleaning the product*).
2. Apply the recommended disinfectant solution. Do not soak.
3. Allow the product to air dry. Follow the instructions and exposure time as directed by the manufacturer of the chemical products.

Using an unsuitable agent, an excessive amount of an agent, an exposure time that is too long or failing to neutralise an agent sufficiently may damage the equipment. All information provided by the manufacturer of the agent used is decisive and must be observed at all times.

7.4 Final treatment

Residues of cleaning agents must be removed. After disinfection, always wipe the product with a clean, water-soaked cloth and allow it to air dry.

Residues of cleaning agents are irritating. They can injure users and damage the product.

8 Maintenance

The operator and his personnel can carry out the following tasks if the instructions are followed. It is the operator's responsibility to employ suitable, qualified and trained personnel carry this out.

8.1 Maintenance

The handle extension is maintenance-free.

Carry out a visual inspection once a week (Section 8.2 *Inspection*).

8.2 Inspection

We recommend cleaning and inspecting the handle extension at the same time. If possible, combine with existing routines (e.g. when cleaning the vehicle).

1. Check the handle extension for any damage as described in Section 8.3.
2. Check the locking button is working properly.
3. Check the system is firmly attached to the clamps on the wheeled stretcher frame.
4. Check that all screws are lying flush.

Inspections must be made weekly.

WARNING

The handle extension must be fully locked. Otherwise, there is a risk of the handle bar detaching from the stretcher, and the user losing control while manoeuvring the stretcher, or the handle bar becoming a projectile in the event of a traffic accident.

8.3 Damage

CAUTION

Damage to the system can create a safety hazard during use. Any damage must be immediately reported to the technical manager.

The handle extension must be withdrawn from use in the case of the following damage:

- Bent handle bar
- Material damage (sharp edges or cracks forming)
- Safety precautions no longer legible
- Missing or loose screws
- Damaged locking mechanism
- Cracks
- Other damage

WARNING

If the wheeled stretcher is withdrawn from use because of an accident, the handle extension attached to it must **not** be reused.

This does not apply to visual defects (e.g. scratches) that cannot cause injury.

8.4 Troubleshooting

Users or technical managers can carry out troubleshooting. Use Table 5: *Troubleshooting*. If the problem cannot be solved or a repair is needed (Section 8.5 *Maintenance and repairs*) please contact us.

Problem	Possible cause	Troubleshooting
Handle bar cannot be inserted	Handlebar stuck in recess	Remove bar and reinsert (Section 6.5 <i>Removing and inserting the handle bar</i>)
	Handle bar bent	Hold the handlebar at both ends and carefully bend back into shape Check from time to time if the handle bar can be reinserted
	Foreign object in the recess	Remove foreign object
Handle extension cannot be locked	Handlebar not fully inserted	Hold the handle bar at both ends and insert it fully If the bar cannot be inserted any further, the bar may be stuck
	Handlebar stuck in recess	See above
	Foreign object in the hole of the handle bar	Remove foreign object
	Locking mechanism defective	Replace locking mechanism

Table 5: *Troubleshooting*

8.5 Maintenance and repairs

The components below can be replaced by the user or operator with original spare parts (Section 12.1 *List of individual parts*). These can be ordered upon request. Repairs are prohibited.

- Handle bar (Item no. fEX201005)
- Quick lock (Item no.: fEX201006)

Only replacement parts made by fEX Engineering GmbH may be used.
Using other parts may be dangerous.

If a locking mechanism needs to be replaced:

1. Unscrew the damaged locking mechanism using a 12 mm wrench.
2. Apply a drop of medium-strength screw locking to the lower threads of the new locking mechanism.
3. Take care to avoid getting any screw locking in the moving parts of the locking mechanism. Wipe away any excess screw locking.
4. Screw the locking mechanism into the outer clamp. Tighten the locking mechanism using a 12 mm torque wrench at 25 Nm.
5. Open the quick-release fastener (red button up; locking pin pulled up). Leave the quick-release fastener in the open position until the screw locking has cured (30 minutes).

The Instructions for use are decisive for all work on the product. Once damaged components have been replaced, the system can be put back into use.

fEX Engineering GmbH refuses to supply spare parts for any mechanisms that have exceeded the end of their service life.

9 Dismantling, recycling, disposal

9.1 Dismantling

The product must be cleaned and disinfected before being dismantled. Components that have been dismantled or damaged may have sharp edges.

The operator and his personnel may dismantle the product. It is the operator's responsibility to have qualified personnel carry this out. The product is dismantled in the reverse order of its assembly.

9.2 Recycling, disposal

When disposing of the product, always follow **current local** recommendations, and/or statutory provisions on environmental protection and risk relating to recycling or the disposal of products at the end of their useful life.

9.3 Disposal through fEX Engineering GmbH

If a product must be withdrawn from use due to damage or other reasons, the handle extension (like all fEX Engineering GmbH products) can be returned to the company for free disposal.

10 Risk management

10.1 Foreseeable misuse

CAUTION

The system must only be used for its intended purpose.

The following foreseeable misuses are prohibited:

- Use without locking or without full locking
- Use of another product than the Stryker *Power-PRO XT*
- Incorrect cleaning or disinfection
- Modification or repair measures by the user that are contrary to these Instructions for use
- Failure to carry out regular inspections
- Use of spare parts made by another manufacturer
- Incorrect installation
- Hanging portable devices on the handle bar
- Using the bar to lift the stretcher, with or without a patient on it
- Using it to sit on
- Other uses contrary to the Instructions for use

Any improper use renders any liability or warranty from the manufacturer null and void.

10.2 Implemented risk control measures

As part of risk analysis, risk control measures to minimise risk have been defined and implemented, in accordance with the Risk Management Plan.

As far as possible, risk control measures are based on compliance with the applicable norms and standards. Technical control measures and monitoring measures have also been taken, and their effectiveness successfully verified.

10.3 Residual risk

Table 10 shows other factors to consider in order to minimise residual risk.

Residual risk	Measures to minimise residual risk
Objects become projectiles inside a vehicle	• Make sure the handle bar is fully locked • Make sure the product is installed correctly • Only use replacement parts from fEX Engineering GmbH
Handle comes off stretcher when manoeuvring	
Sharp edges	• Make sure the product is handled and used correctly • Carry out weekly inspections • Withdraw damaged components from use immediately

Table 6: Measures to minimise residual risk

11 Performance data

11.1 Product conformity

The product meets the following regulatory requirements:

- DIN EN 1789:2020-12 Medical vehicles and their equipment - Road ambulances
- (EU) 2017/745 (MDR) REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on Medical devices
- ISO 14971:2022-04 Medical devices - Application of risk management to medical devices
- IEC 62366-1:2021-08 Medical devices – Part 1: Application of usability engineering to medical devices
- DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes

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