International Instructions for Use **fEX Equipment Rail**

Languages:

EN, DA, DE,
ES, FI, FR, HU,
IT, NL, NO, SL,
SV

fEX Engineering GmbH

Höltenweg 49
48155 Münster
Germany

Fourth edition

Publication date 22/07/2024

Foreword

Thank you for purchasing an **fEX equipment rail** from fEX Engineering GmbH.

The fEX equipment rail is an extension for medical products such as wheeled stretchers, patient couches, bed frames and rescue vehicles, according to standard DIN EN ISO 1789. This rail can be used to attach, support and move additional medical devices to the products.

WARNING

Improper handling or any unintended use may cause incur risk or cause damage. Please therefore read and carefully follow these Instructions for use. Retain them for reference at all times. To prevent personal injury or damage to property, please observe all safety precautions.

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster
Germany
engineering@fex.group
+49 (0)251 508 537 -13

Contents

1	How to use these Instructions for use.....	1
1.1	Feedback on these Instructions for use	1
1.2	Explanation of the warning notices used in these Instructions for use	1
1.3	Abbreviations	2
2	Intended purpose	3
2.1	Intended use.....	3
2.2	Intended medical purpose	3
2.3	Clinical benefits	3
2.4	Indications	3
2.5	Contraindications	3
2.6	Intended target patient group	3
2.7	Intended use with other products.....	3
2.8	Intended users.....	4
2.9	Intended environment for use	4
3	Important safety notice	5
3.1	What to do in the event of damage	5
3.2	Safety notice on load capacity.....	5
3.3	Reporting obligation.....	6
4	Product description	7
4.1	Product overview	7
4.2	Symbols used.....	8
4.3	How it works.....	8
4.4	Use in combination with other products	9
4.5	Delivery.....	10
4.6	Modifications.....	10
4.7	Requirements for storage,handling and environmental conditions	10
4.8	Product service life	10
4.9	Utility model/patent.....	10
5	Assembly and putting into service	11
5.1	Requirements for personnelinvolved.....	11
5.2	General	11
5.3	Stryker Prime Series - assembly with handle clamp.....	12
5.4	Stryker Prime Series – assembly with arm	13
5.5	Stryker Prime X Series – assembly with side clamp	14
5.6	Stryker Fluoroscopy Stretcher	15

5.7	Stryker SV2 hospital bed	16
5.8	Stryker ProCuity hospital bed	17
5.9	Stryker M1 Roll-In System wheeled stretcher.....	18
5.10	Stryker Power-PRO XT wheeled stretcher.....	19
5.11	Stryker Power-PRO 2 wheeled stretcher.....	20
5.12	Wall bracket.....	21
5.13	To check that all variants have been assembled or installed correctly.....	21
6	Using the product	22
6.1	Operator's duty of care	22
6.2	User's duty of care.....	22
6.3	Instructions for safe operation.....	22
7	Reprocessing: cleaning, disinfection	23
7.1	Cleaning the product	23
7.2	Disinfectants	23
7.3	Disinfecting the product:.....	24
7.4	Final treatment.....	24
8	Maintenance.....	25
8.1	Maintenance.....	25
8.2	Inspection	25
8.3	Damage.....	25
8.4	Maintenance and repairs	26
9	Dismantling, recycling, disposal	27
9.1	Dismantling.....	27
9.2	Recycling, disposal.....	27
9.3	Disposal through fEX Engineering GmbH	27
10	Risk management	28
10.1	Foreseeable misuse	28
10.2	Implemented risk control measures	28
10.3	Residual risk.....	29
11	Performance data.....	30
11.1	Product conformity	30
12	Technical data.....	31
12.1	List of individual parts	31

1 How to use these Instructions for use

If you have any questions about the contents of these Instructions for use or how to use the product, please contact us. You can request a paper copy of these Instructions for use or download the current, electronic version from our website:

<https://www.fex.group/downloads>

Should you need further assistance after carefully reading these Instructions for use and additional information, please contact your authorised distributor or contact us directly.

1.1 Feedback on these Instructions for use

Your opinion is important to us. Please let us know if there is anything you would change about these Instructions for use (engineering@fex.group). We will analyse your feedback and take it into account in the next version of these Instructions for use.

1.2 Explanation of the warning notices used in these Instructions for use

Signal word	Description
 WARNING	Potentially hazardous situation resulting in severe injury or death.
 CAUTION	Potential hazard resulting in minor or moderate injury.
NOTE	User error can damage the product.

Table 1: Explanation of the warning notices used in these Instructions for use

1.3 Abbreviations

Abbreviation	Definition
Item no.	Item number
BfArM	Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte)
BHCS	Button head cap screw
HFB	Hex flange bolt
MDR	Medical Device Regulation
MPSV	Medical Devices Safety Plan Ordinance (Medizinprodukte-Sicherheitsplanverordnung)
N, Nm	Newton, Newton-meter
PL	Patient Left
PLF / PRF	Patient Left Foot-end / Patient Right Foot-end
PLH / PRH	Patient left head-end / Patient right head-end
PR	Patient right
SHCS	Socket head cap screw
SSSFP	Socket set screw flat point
UDI	Unique Device Identification
UDI-DI	Unique Device Identification – Device Identifier
UDI-PI	Unique Device Identification – Production Identifier
UNC	Unified national coarse thread

Table 2: Abbreviations

2 Intended purpose

2.1 Intended use

The fEX equipment rail has been designed as an extension for medical products according to standard DIN EN ISO 19054. It can be attached to stretchers, beds and in rescue vehicles according to standard DIN EN ISO 1789, and used to support and transport medical devices.

2.2 Intended medical purpose

The fEX equipment rail is a Class I, mobile, non-invasive, non-active medical device.

2.3 Clinical benefits

The fEX equipment rail is classed as an accessory and has no inherent clinical benefit. It is designed to allow the intended purpose of any medical devices attached to the rail to be achieved.

2.4 Indications

No indications have been specified for the fEX equipment rail.

2.5 Contraindications

There are no known contraindications for the fEX equipment rail.

2.6 Intended target patient group

The patient population is actually using the medical device attached to the rail.

No direct contact with the patient is anticipated under normal circumstances.

2.7 Intended use with other products

fEX equipment rails are extensions for medical products, and should only be mounted on and used with products approved by the manufacturer. Approved products are listed in Section 4.4 *Use in combination with other products*.

2.8 Intended users

fEX equipment rails are intended for use in a professional healthcare environment. The users, and the qualifications they require, are specified accordingly.

User groups	Nursing staff, emergency services personnel, doctors (primary users)	Medical technicians, technical managers, cleaning staff (secondary users)
Training, certification	Medically trained staff, e.g. emergency paramedics, emergency medical technicians, rescue workers, emergency services support staff, doctors, or qualified nurses	Medical technicians (training/academic qualifications) or equivalent, or accreditation as a cleaner in a professional healthcare environment

Table 3: Intended users

2.9 Intended environment for use

fEX equipment rails are intended for use by primary users in the following environments:

- Use in proximity to patients and therefore in a professional healthcare environment
- Use in a mobile environment, e.g. in emergency services vehicles
- Use outdoors for rescue operations

fEX equipment rails do not require any specific physical settings.

3 Important safety notice

NOTE

Please read and observe the instructions in these Instructions for use to use this system safely.

CAUTION

Only use battery-powered medical devices on *M1 Roll-In*, *Power-PRO XT* and *Power-PRO 2* wheeled stretchers.

If a wheeled stretcher can no longer be used because of an accident or similar incident, the same applies to the fEX equipment rail fitted to it.

WARNING

Only use rail clamps, rail brackets and devices permitted in accordance with standard DIN EN ISO 19054. Otherwise, there is a danger of parts flying off.

3.1 What to do in the event of damage

fEX extension rails must be withdrawn from service if they are damaged (Section 8.3 *Damage*). Repairing damaged components is prohibited – these must be replaced (Section 8.4 *Maintenance and repairs*).

Visible damage must be immediately reported to the technical manager.

3.2 Safety notice on load capacity

The maximum permitted load is engraved on the rail and must not be exceeded:

- Wheeled stretchers (e.g. *Power-PRO 2*): 15 kg
- Other products used with the rail: 20 kg

WARNING

No devices may be fitted to the wheeled stretcher during a journey in a vehicle. These must be stored or secured elsewhere (e.g. wall bracket).

3.3 Reporting obligation

All serious incidents which occur in connection with the product must be reported to the manufacturer and the competent authority in the EU Member State where you are established.

In Germany, the competent authority is BfArM. Up-to-date contact details can be found on the BfArM website: <https://www.bfarm.de>

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster
Germany
+49 (0)251 508 537 -13
engineering@fex.group

4 Product description

4.1 Product overview

The rail is made of stainless steel, and supported by two aluminium rail brackets. fEX equipment rails feature a standard design and therefore require a suitable mounting accessory.

Figure 1: Product photo

The accessory required depends on the product being attached to the equipment rail. Mounting accessories must only be used for approved products.

4.2 Symbols used


Figure 2: Notices on the product

Notice	Meaning
Max. 15 kg or Max. 20 kg	Maximum load
fEX GROUP MEDICAL SERVICE & ENGINEERING	Manufacturer's logo
MD	Medical device
CE	Conformity marking
	Read Instructions for use
	Manufacturer, address and contact details
	UDI – Unique Device Identification: (00) UDI-DI (Device Identifier; (10) Batch no.; (11) Date of manufacture

Table 4: Notices on the product

4.3 How it works







The fEX equipment rail has been designed in accordance with standard DIN EN ISO 19054 3.9. Medical devices can be attached to it and supported by it.

The rail consists of a hollow profile held by two rail brackets. An additional mounting accessory is required to mount the rail on stretchers and beds, and in emergency services vehicles (Section 4.4 *Use in combination with other products*).

The addition of rail clamps, in accordance with standard DIN EN ISO 19054, Subsection 3.10, means that the fEX equipment rails form a rail system in accordance with standard DIN EN ISO 19054 3.15.

4.4 Use in combination with other products

fEX equipment rails are intended for use in combination with other products. The rail is approved for the products listed in Table 5: *Product variants* and can be used with the mounting accessories listed there.

Variant	Item number	Connector
V1: Stryker Prime Series with handle clamp	fEX101026 (26" stretcher) fEX101030 (30" stretcher)	
V2: Stryker Prime Series with arm	fEX102026 (26" stretcher) fEX102030 (30" stretcher)	
V3: Stryker Prime X Series with side clamp	fEX103001 (one-sided, short rail) fEX103002 (one-sided, long rail) fEX103003 (two-sided, short rail) fEX103004 (two-sided, long rail)	
V4: Stryker Fluoroscopy with arm	fEX104001	
V5: Stryker SV2	fEX105001	
V6: Stryker ProCuity	fEX106001	
V7: Stryker M1-Roll-in-System	fEX1070PL (patient left, foot-end) fEX1070PR (patient right, foot-end)	




V8: Stryker Power-PRO XT	fEX1080PL (patient left, foot-end) fEX1080PR (patient right, foot-end)	
V9: Stryker Power-PRO 2	fEX1090PL (patient left, foot-end) fEX1090PR (patient right, foot-end)	
V10: Wall bracket	fEX110003 (clamping range 300 mm) fEX110004 (clamping range 400 mm) fEX110005 (clamping range 500 mm) fEX110006 (clamping range 600 mm)	

Table 5: Product variants

4.5 Delivery

fEX equipment rails are delivered as a complete assembly, in non-sterile, non-disinfected condition.

An electronic copy of the Instructions for use can be found via the following link:

<https://www.fex.group/downloads/>

A sticker on the product shows a QR code and the above link, both of which will take you to the electronic version of the Instructions for use.

4.6 Modifications

Unauthorised modifications to this product are prohibited.

4.7 Requirements for storage, handling and environmental conditions

The product requires no specific storage or handling conditions.

4.8 Product service life

Under normal usage conditions and with proper maintenance, the service life is **7 years**.

The date for the product's service life is based on the delivery date from fEX Engineering.

4.9 Utility model/patent

Registered utility model: No. 20 2021 105 211

5 Assembly and putting into service

5.1 Requirements for personnel involved

The operator and his personnel can assemble the product by following the instructions below. It is the operator's responsibility to have qualified personnel carry this out.

5.2 General

fEX equipment rails are delivered as a complete assembly as per the order (for the product to be used with the rail, e.g. Stryker *Power-PRO XT*), i.e. the package will include rails, brackets and mounting accessories. To ensure traceability, in accordance with the MDR, this assembly must be installed as a coherent unit as it is delivered.

A torque wrench must be used for assembly; the specified torques must be observed.

Medium-strength screw locking is strongly recommended. The required curing time depends on the information given by the screw locking manufacturer.

fEX equipment rails are delivered in non-disinfected condition. You must disinfect the handle extension before first use. For disinfection instructions, see Section 7 *Reprocessing: cleaning and disinfection*.

The following general information applies to all variants:

- Rails are generally delivered with the rail brackets pre-assembled
- If necessary, slide the standard rail into the rail brackets
- Markings on the standard rail must be visible to the user, but not to the patient (Figure 2: *Rail with rail bracket*).


Figure 3: Rail with rail brackets

5.3 Stryker Prime Series - assembly with handle clamp

This is how to assemble the rail system on patient stretchers in the Stryker *Prime Series* with clamps on the handles. Secure clamping can only be guaranteed for handles with a diameter of 32 mm to 34 mm.


Figure 4: Assembly instructions for Prime handle clamp


Figure 5: Prime handle

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Using screws (2) UNC ¼ -20 x ½ (SHCS), screw the *rail bracket* (1) to the *front clamp* (3) using a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 as required for the other side.
5. Hold the rail system by the handle connector. To avoid blocking the folding mechanism, the clamp in the lower area must be in position (7).
6. Screw on the *rear clamp* (5) using screws (6) UNC ¼ -20 x 1 ¼ (SHCS) and a tightening torque of 12 Nm.

5.4 Stryker Prime Series – assembly with arm

This is how to assemble the rail system on patient stretchers in the Stryker *Prime Series* without handles and with the *arm* accessory. Use the free moving sleeve (6) here.


Figure 6: Assembly instructions for Prime arm


Figure 7: Prime connector

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) onto the *arm* (3) using screws (2) UNC ¼ -20 x ¾ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Insert the shaft of the *arm* (3) into the sleeve (6) of the support.
6. Screw the *flange screw* (5) UNC ¼ -20 x ½ (HFB) into the shaft of the *arm* underneath with a tightening torque of 12 Nm.

5.5 Stryker Prime X Series – assembly with side clamp

This is how to assemble the rail system on patient stretchers in the *Prime-X Series* made by *Stryker* by clamping it to the side of the bar around the lying surface of the stretcher.


Figure 8: Assembly instructions for Prime X


Figure 9: Prime X Rohr

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw on the *rail bracket* (1) with screws (2) UNC ¼ -20 x ½ (SHCS) and the *outside clamp* (3) with a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Hold the rail system on the bar (7) around the lying surface in position (8/9) inside the handle recess.
6. Screw the *inside clamp* (5) to the *outside clamp* using screws (6) UNC ¼ -20 x 1 ¼ (SHCS) with a tightening torque of 12 Nm. Make sure that the flat side of the *inside clamp* (5) is facing up towards the lying surface.

5.6 Stryker Fluoroscopy Stretcher

This is how to fit the equipment rail to the free moving connector on the *Fluoroscopy Stretcher* made by *Stryker* when it is already fitted with the *arm* accessory.


Figure 10: Assembly instructions for Fluoroscopy


Figure 11: Fluoroscopy connector

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) to the *arm* (3) using screws (2) UNC ¼ -20 x ¾ (SHCS) using a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Insert the shaft of the *arm* (3) into the connectors (6) for the support.
6. Screw in flange screws (5) UNC ¼ -20 x ½ (HFB) with a tightening torque of 12 Nm under the shaft of the *arm*.

5.7 Stryker SV2 hospital bed

This is how to fit the foot-end equipment rail on the SV2 patient bed made by *Stryker*.


Figure 12: Assembly instructions for SV2


Figure 13: SV2 screw connections

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) and the *holding plate* (3) together using screws (2) UNC $\frac{1}{4}$ -20 x $\frac{3}{4}$ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Remove the panel from the foot-end of the bed and remove the four original *Stryker* screws and nuts (9).
6. Place the rail system over the screw holes left empty by step 5.
7. Insert screws (5) UNC #10 -20 x $\frac{3}{4}$ through the *holding plate* (3) and through the *spacer* (8).
8. Put on the washers (7) and screw in the nuts (6) with a tightening torque of 8 Nm.
9. Finally, replace the *Stryker* panel.

5.8 Stryker ProCuity hospital bed

This is how to fit the foot-end equipment rail on the *ProCuity* patient bed made by *Stryker*.


Figure 14: Assembly instructions for the ProCuity


Figure 15: ProCuity (foot-end)

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) onto the *arm* (3) using screws (2) UNC ¼ -20 x ½ with a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Remove the panel from the foot-end of the bed (6).
6. Insert the rail system into the connectors (7/8).
7. Screw in the safety screws (5) UNC ¼ -20 x 5/16 with a tightening torque of 7 Nm.
8. Finally, replace the *Stryker* panel.

5.9 Stryker M1 Roll-In System wheeled stretcher

This is how to fit the equipment rail to an *M1 Roll-In System* rescue stretcher made by Stryker.


Figure 16: Assembly instructions for the M1 Roll-In


Figure 17: M1 Roll-In IV Pole and Clip

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) onto the *holding plate* (3) using screws (2) UNC ¼ -20 x ¾ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Unscrew the four original *Stryker* screws (5) and remove the *IV Pole* (6) and its retaining clip (7) (if used).
6. Place the rail system over the screw holes left empty by step 5.
7. If the *IV Pole* (6) and retaining clip (7) are used, put these back **on top** of the *holding plate* (3).
8. Screw in the original *Stryker* screws (5) UNC ¼ -20 x 1 (BHCS) using a tightening torque of 12 Nm.

5.10 Stryker Power-PRO XT wheeled stretcher

This is how to mount the rail system on the side of the Power-PRO XT rescue stretcher made by Stryker using the *XT holding plate* accessory.


Figure 18: Assembly instructions for the Power-PRO XT


Figure 19: Power-PRO XT IV Pole and Clip

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) and the *holding plate* (3) together using screws (2) UNC ¼ -20 x ¾ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 for the other side.
5. Unscrew the four original *Stryker* screws (5) and remove the *IV Pole* (6) and the retaining clip (7) (if used).
6. Place the rail system over the screw holes left empty by step 5.
7. If the *IV Pole* (6) and retaining clip (7) are used, put these back **on top** of the *holding plate* (3).
8. Screw in the original *Stryker* screws (5) UNC ¼ -20 x 1 (BHCS) using a tightening torque of 12 Nm.

5.11 Stryker Power-PRO 2 wheeled stretcher

This is how to attach the rail system to the side of the *Power-PRO 2* rescue stretcher made by *Stryker* with the *holding plate* accessory.


Figure 20: Assembly instructions for the Power-PRO 2


Figure 21: Power-PRO 2 IV Pole and Clip

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) onto the holding plate (3 or 6) using screws (2) UNC ¼-20x ½ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 as required for the other side.
5. To use **with** the *IV Pole* (7):
 - a. Unscrew the six original *Stryker* screws (9) and remove the *IV Pole* (7).
 - b. Place the rail system **on top** of the *retaining clip* (8).
 - c. Place the *IV Pole* (7) **on top** of the *holding plate* (6) for the rail system.
 - d. Screw in the accompanying screws (5) M6x16 (countersunk head screw) with a tightening torque of 12 Nm.
6. To use **without** the *IV Pole* (7):
 - a. Place the rail system over the empty screw holes.
 - b. Screw in the accompanying screws (5) M6x16 (countersunk head screw) with a tightening torque of 12 Nm.

5.12 Wall bracket

This is how to mount the rail system on the wall of a building or in a rescue vehicle, according to DIN EN 1789, and using the *flange* accessory. M8 threaded holes or riveting nuts must be used to attach the wall bracket.


Figure 22: Assembly instructions for the wall bracket

Follow the steps below:

1. Apply a drop of medium-strength screw locking to the first (lowest) threads on all screws.
2. Screw the *rail bracket* (1) and the *flange* (3) together using screws (2) UNC ¼-20x ¾ (SHCS) and a tightening torque of 12 Nm.
3. Insert the *rail* (4) into the *rail bracket* (1).
4. Repeat steps 1 and 2 as required for the other side.
5. Hold the system against the wall or in front of the threaded holes or riveting nuts.
6. Screw in the accompanying screws (5) M8x25 (countersunk head screw) with a tightening torque of 25 Nm.

5.13 To check that all variants have been assembled or installed correctly

WARNING

Use the check list below to make sure variants have been assembled or installed correctly. If any of these points do not apply, the product must be reassembled or reinstalled.

- The maximum load engraving should be visible to the user in the direction of the patient
- The rail should sit inside the rail holder with no play
- All necessary screws must have been used (no through holes are visible)
- All screws should be flush and tightly screwed (use a torque wrench for assembly).
- The accessories used cannot be loosened/moved

Once assembled or installed properly, the rail system is ready for immediate use

6 Using the product

6.1 Operator's duty of care

To ensure safe operation, the operator must meet the following obligations:

- Personnel must be trained to use the rail system installed
- The Instructions for use must be read and followed
- Products must be assembled or installed properly
- Only suitable, qualified, trained specialist personnel may use the product
- All serious incidents (damage, injuries, etc.) must be reported to the manufacturer and the competent authorities

6.2 User's duty of care

To ensure safe operation, the user must meet the following obligations:

- The Instructions for use must be read and followed
- The product must be used as intended
- Damage must be reported to the technical manager

6.3 Instructions for safe operation

- Do not exceed the maximum load
- No devices may be fitted to the rail of the wheeled stretcher during a journey in a vehicle
- Only use the product for its intended purpose
- Do not use to pull, push or lift
- Only use rail clamps, rail holders and devices with a holder that is in accordance with standard DIN EN ISO 19054
- Only use battery-powered medical devices on the *M1 Roll-in-System*, *Power-PRO XT* and *Power-PRO 2* wheeled stretchers made by Stryker.

7 Reprocessing: cleaning, disinfection

The fEX equipment rail is a reusable medical device. Cleaning and disinfection must follow the process in the Instructions for use, and be carried out after every potential contamination, or at least once a month.

We recommend cleaning and inspecting the fEX equipment rail at the same time. If possible, combine with existing weekly routines (e.g. when cleaning the vehicle).

7.1 Cleaning the product

The product can be cleaned using high pressure.

1. Wipe the product's surfaces with a solution of mild soap and water.
2. Dry all product surfaces thoroughly.

7.2 Disinfectants

Use the following agents for disinfection:

- Quaternary cleaners containing less than 3% glycol ether
- Chlorinated bleaching agents (5.25 % – less than 1 part bleaching agent to 100 parts of water)
- Phenolic cleaners (active ingredient: o-phenylphenol)
- ≤ 21% isopropanol

Disinfectants made from quaternary ammonium compounds must be used in accordance with the recommended instructions for use.

WARNING

Disinfectants containing any of the following basic ingredients are not approved as they can damage the system:

- Hydrogen peroxide
- (Per-/Chlor-) acetic acid

7.3 Disinfecting the product:

1. Remove any visible contaminants from the product (Section 7.1 *Cleaning the product*).
2. Apply the recommended disinfectant solution. Do not soak.
3. Allow the product to air dry. Follow the instructions and exposure time as directed by the manufacturer of the chemical products.

Using an unsuitable agent, an excessive amount of an agent, an exposure times that is too long or failing to neutralise an agent sufficiently may damage the equipment. All information provided by the manufacturer of the agent used is decisive and must be observed at all times.

7.4 Final treatment

Residues of cleaning agents must be removed. After disinfection, always wipe the product with a clean, water-soaked cloth and allow it to air dry.

Residues of cleaning agents are corrosive and can injure users and damage the product.

8 Maintenance

The operator and his personnel can carry out the following tasks if the instructions are followed. It is the operator's responsibility to employ suitable, qualified and trained personnel carry this out.

8.1 Maintenance

The rail system is maintenance-free.

Carry out a visual inspection once a week (Section 8.2 *Inspection*).

8.2 Inspection

We recommend inspecting the fEX equipment rail when cleaning it. If possible, combine with existing routines (e.g. when cleaning the vehicle).

1. Check the fEX equipment rail for any damage as described in Section 8.3.
2. Check the system is firmly attached to the clamps on the wheeled stretcher frame.
3. Check for any screws that are not lying flush.

Inspections should be made weekly.

8.3 Damage

CAUTION

Damage to the system can create a safety hazard during use. Any damage must be immediately reported to the technical manager.

fEX equipment rails must be withdrawn from use in the case of the following damage:

- Bent plates/components
- Bent equipment rails
- Material damage (sharp edges or cracks forming)
- Safety precautions no longer legible
- Missing or loose screws
- Other damage

WARNING

If the wheel stretcher is withdrawn from use because of an accident, the equipment rails attached to it must **not** be reused.

This does not apply to visual defects such as superficial scratches that cannot cause injury.

8.4 Maintenance and repairs

All components can be replaced by the user or operator with original spare parts. The components for all variants and their item numbers are listed in Section 12.1 *List of individual parts*. These can be ordered upon request.

WARNING

Only replacement parts made by fEX Engineering GmbH may be used.
Using other parts may be dangerous.
Repairs are prohibited.

The Instructions for use are decisive for all work on the product. Once damaged components have been replaced, the system can be put back into use.

fEX Engineering GmbH refuses to supply spare parts for any products that have reached the end of their service life.

9 Dismantling, recycling, disposal

9.1 Dismantling

**CAUTION**

The product must be cleaned and disinfected before being dismantled.
Components that have been dismantled or damaged may have sharp edges.

The operator and his personnel may dismantle the product. It is the operator's responsibility to have qualified personnel carry this out. The product is dismantled in the reverse order of its assembly.

9.2 Recycling, disposal

When disposing of the product, always follow **current local** recommendations, and/or statutory provisions on environmental protection and risk relating to recycling or the disposal of devices at the end of their useful life.

9.3 Disposal through fEX Engineering GmbH

If a product must be withdrawn from use due to damage or other reasons, the equipment rail (like all fEX Engineering GmbH products) can be returned to the company for free disposal.

10 Risk management

10.1 Foreseeable misuse

CAUTION

The system must only be used for its intended purpose.

The following foreseeable misuses are prohibited:

- Exceeding the maximum load
- Use as a handle for lifting, pushing or pulling
- Use as a seat or for leaning
- Modification or repair measures by the user that are contrary to these Instructions for use
- Incorrect cleaning or disinfection
- Failure to carry out regular inspections
- Incorrect installation
- Attaching devices without a holder that is suitable for equipment rails (according to standard DIN EN ISO 19054) e.g. a rail clamp
- Other uses contrary to the Instructions for use

Any improper use renders any liability or warranty from the manufacturer null and void.

10.2 Implemented risk control measures

As part of risk analysis, risk control measures to minimise risk have been defined and implemented, in accordance with the Risk Management Plan.

As far as possible, risk control measures are based on compliance with the applicable norms and standards. Technical control measures and monitoring measures have also been taken, and their effectiveness successfully verified.

10.3 Residual risk

Table 6 shows other factors to consider in order to minimise residual risk.

Residual risk	Measures to minimise residual risk
Exceeding the maximum load	• Ensuring users are trained • Checking legibility of instructions, otherwise replacing equipment rails
Objects become projectiles inside a vehicle	• Only use rail clamps, rail holders and devices with a holder that are in accordance with standard DIN EN ISO 19054. • Make sure the product is installed correctly • Make sure the product is handled and used correctly • Carry out weekly inspections • Withdraw damaged components from use immediately • Only use replacement parts from fEX Engineering GmbH
Sharp edges	• Make sure the product is handled and used correctly • Carry out weekly inspections • Withdraw damaged components from use immediately
Devices fall off the rail and injure users	• Observe the maximum load • Only use rail clamps rail brackets and devices permitted in accordance with standard DIN EN ISO 19054
Potential electrical charge	• Make sure to only use battery-powered medical devices on <i>M1 Roll-In</i>, <i>Power-PRO XT</i> and <i>Power-PRO 2</i> wheeled stretchers. • Ensure that the insulation on attached medical devices is intact. If this insulation fails and the devices are not earthed through the frame (of the stretcher, bed, etc.) the rail may be live.
Potential transmission of pathogens	• Ensure the product is cleaned and disinfected properly according to the procedure in the Instructions for use. • Clean and disinfect fEX equipment rails after every possible contamination.

Table 6: Measures to minimise residual risk

11 Performance data

11.1 Product conformity

The product meets the following regulatory requirements:

- DIN EN ISO 19054:2006 + A1:2016 Rail systems for supporting medical equipment
- DIN EN 1789:2020-12 Medical vehicles and their equipment - Road ambulances
- (EU) 2017/745 (MDR) REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on Medical devices
- ISO 14971:2022-04 Medical devices - Application of risk management to medical devices
- IEC 62366-1:2021-08 Medical devices – Part 1: Application of usability engineering to medical devices
- DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes

12 Technical data

12.1 List of individual parts

V.	Assembly	Components	Quantity per set	Characteristic / designation	Value
	Equipment rail		1	Dimensions (height x depth)	25 mm x 10 mm
				Weight	Depending on variant (0.8 g/mm)
				Material	Stainless steel (brushed)
	Rail bracket (Item no.: 100003)		2	Dimensions (width x height x depth)	62 mm x 31 mm x 25 mm
				Weight	79 g
				Material	Aluminium (anodised)
1	Prime Series handle clamp Item no.: fEX101026/ fEX101030	Equipment rail (Item no.: 101003 / 101004)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 587 mm (26") 25 mm x 10 mm x 688 mm (30")
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000SHCS03)	4	SHCS	UNC ¼-20x 1 ¼
		Equipment rail bracket and clamps (Item no.: fEX101001)	2	Dimensions (width x height x depth)	62 x 31 x 75 mm
				Weight	190 g
				Material	Aluminium (anodised)
2	Prime Series arm Item no.: fEX102026/ fEX102030	Equipment rail (Item no.: 102003 / 102004)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 587 mm (26") 25 mm x 10 mm x 688 mm (30")
		Screws (Item no.: 000SHCS02)	4	SHCS	UNC ¼-20x ¾
		Screws (Item no.: 000HFB001)	2	HFB	UNC ¼-20x ½
		Equipment rail bracket and PL arm (Item no.: fEX1020PL)	1	Dimensions (width x height x depth)	91 mm x 128 mm x 158 mm
				Weight	325 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel
		Equipment rail bracket and PR arm (Item no.: fEX1020PR)	1	Dimensions (width x height x depth)	91 mm x 128 mm x 158 mm
				Weight	325 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel
3	Prime X Series side bar clamp Item no.: fEX103001/ fEX103002/ fEX103003/ fEX103004	Equipment rail (Item no.: 103005 / 103006)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 285 mm 25 mm x 10 mm x 745 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000SHCS03)	8	SHCS	UNC ¼-20x 1 ¼
		Equipment rail bracket and outer lock/inner bar (Item no.: fEX103HMS)	1	Dimensions (width x height x depth)	62 mm x 47 mm x 108 mm
				Weight	260 g
				Material	Aluminium (anodised)
		Equipment rail bracket and inner/outer bar clamp (Item no.: fEX103F00)	1	Dimensions (width x height x depth)	62 mm x 47 mm x 108 mm
				Weight	265 g
				Material	Aluminium (anodised)
		Equipment rail bracket and inner/outer bar clamp (Item no.: fEX103H00)	2	Dimensions (width x height x depth)	62 mm x 47 mm x 108 mm
				Weight	265 g
				Material	Aluminium (anodised)

V.	Assembly	Components	Quantity per set	Characteristic / designation	Value
4	Fluoroscopy arm Item no.: fEX104001	Equipment rail (Item no.: 104003)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 645 mm
		Screws (Item no.: 000SHCS02)	4	SHCS	UNC ¼-20x ¾
		Screws (Item no.: 000HFB001)	2	HFB	UNC ¼-20x ½
		Equipment rail bracket and PL arm (Item no.: fEX1040PL)	1	Dimensions (width x height x depth)	91 mm x 78 mm x 158 mm
				Weight	295 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel
		Equipment rail bracket and PR arm (Item no.: fEX1040PR)	1	Dimensions (width x height x depth)	91 mm x 78 mm x 158 mm
				Weight	295 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel
5	SV2 Item no.: fEX105001	Equipment rail (Item no.: 105003)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 711 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000BHCS01)	4	BHCS	UNC #10-20x ¾
		Nut (Item no.: 000TNLN01)	4	TNLN	UNC ¼
		Washer (Item no.: 000FW0001)	4	Washer	¼ inch
		Equipment rail bracket and PR holding plate (Item no.: fEX1050PR)	1	Dimensions (width x height x depth)	62 mm x 43 mm x 111 mm
				Weight	335 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PL holding plate (Item no.: fEX1050PL)	1	Dimensions (width x height x depth)	62 mm x 43 mm x 111 mm
				Weight	335 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Spacer (Item no.: 105002)	2	Dimensions (width x height x depth)	56 mm x 22 mm x 8 mm
				Weight	13 g
				Material	POM
6	ProCuity Item no.: fEX106001	Equipment rail (Item no.: 106003)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 844 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000SSFP1)	4	SSFP	UNC ¼-20x ¾
		Equipment rail bracket and PL arm (Item no.: fEX1060PL)	1	Dimensions (width x height x depth)	62 mm x 31 mm x 150 mm
				Weight	230 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel
		Equipment rail bracket and PR arm (Item no.: fEX1060PR)	1	Dimensions (width x height x depth)	62 mm x 31 mm x 150 mm
				Weight	230 g
				Material	Aluminium (powder coated), aluminium (anodised), stainless steel

V.	Assembly	Components	Quantity per set	Characteristic / designation	Value
7	M1 Roll-In System patient left Item no.: fEX1070PL	Equipment rail (Item no.: 107005)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 378 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼ -20 x ½
		Equipment rail bracket and PLF holding plate (Item no.: fEX107PLF)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PLH holding plate (Item no.: fEX107PLH)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
	M1 Roll-In System patient right Item no.: fEX1070PR	Equipment rail (Item no.: 107005)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 378 mm
		Screws (Item no.: 000SHCS01)		SHCS	UNC ¼ -20 x ½
		Equipment rail bracket and PRF holding plate (Item no.: fEX107PRF)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PRH holding plate (Item no.: fEX107PRH)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
8	Power-PRO XT patient left Item no.: fEX1080PL	Equipment rail (Item no.: 108005)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 333 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Equipment rail bracket and PLF holding plate (Item no.: fEX108PLF)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PLH holding plate (Item no.: fEX108PLH)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
	Power-PRO XT patient right Item no.: fEX1080PR	Equipment rail (Item no.: 108005)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 333 mm
		Screws (Item no.: 000SHCS01)	4	SCHS	UNC ¼-20x ½
		Equipment rail bracket and PRF holding plate (Item no.: fEX108PRF)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PRH holding plate (Item no.: fEX108PRH)	1	Dimensions (width x height x depth)	89 mm x 52 mm x 85 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)

V.	Assembly	Components	Quantity per set	Characteristic / designation	Value
9	Power-PRO 2 patient left Item no.: fEX1090PL	Equipment rail (Item no.: 109001)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 378 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000LKS001)	6	Countersunk head screw	M6 x 16
		Equipment rail bracket and PLH holding plate (Item no.: fEX109PLH)	1	Dimensions (width x height x depth)	83 mm x 52 mm x 86 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PLF holding plate (Item no.: fEX109PLF)	1	Dimensions (width x height x depth)	112 mm x 52 mm x 98 mm
				Weight	275 g
				Material	Mild steel (powder coated), aluminium (anodised)
	Power-PRO 2 patient right Item no.: fEX1090PR	Equipment rail (Item no.: 109001)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 378 mm
		Screws (Item no.: 000SHCS01)	4	SHCS	UNC ¼-20x ½
		Screws (Item no.: 000LKS001)	6	Countersunk head screw	M6 x 16
		Equipment rail bracket and PRH holding plate (Item no.: fEX109PRH)	1	Dimensions (width x height x depth)	83 mm x 52 mm x 86 mm
				Weight	260 g
				Material	Mild steel (powder coated), aluminium (anodised)
		Equipment rail bracket and PRF holding plate (Item no.: fEX109PRF)	1	Dimensions (width x height x depth)	112 mm x 52 mm x 98 mm
				Weight	275 g
				Material	Mild steel (powder coated), aluminium (anodised)
10	Wall assembly Item no.: fEX110003/ fEX110004/ fEX110005/ fEX110006	Equipment rail (Item no.: 110003)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 417 mm
		Equipment rail (Item no.: 110004)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 517 mm
		Equipment rail (Item no.: 110005)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 617 mm
		Equipment rail (Item no.: 110006)	1	Dimensions (width x height x depth)	25 mm x 10 mm x 717 mm
		Screws (Item no.: 000SHCS02)	4	SHCS	UNC ¼-20x ¾
		Screws (Item no.: 000SKS001)	4	Countersunk head screw	M8 x 25
		Equipment rail bracket and flange (Item no.: fEX110002)	2	Dimensions (width x height x depth)	62 mm x 75 mm x 45 mm
				Weight	185 g
				Material	Aluminium (anodised)

Table 7: Device dimensions

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